

WELCH ALLYN ELITE™ RETINOSCOPE



Superior Illumination in a Prestigious New Retinoscope

The Welch Allyn Elite Retinoscope's enhanced light output makes refractions easier than ever by producing one of the brightest, sharpest, and cleanest streaks available. Features such as a convenient external focusing sleeve, sealed optics, and an optional polarizing filter make the Welch Allyn Elite Retinoscope one of the most user-friendly and innovative instruments offered.

WelchAllyn®


Advancing Frontline Care™

SUPERIOR ILLUMINATION FOR UNPARALLELED PERFORMANCE

For more than 85 years, Welch Allyn diagnostic instruments have set the standards for superior performance and value. And with the introduction of the Welch Allyn Elite™ Retinoscope, we've raised the bar once again. This prestigious instrument offers a combination of outstanding features never before seen in a retinoscope. Which means you'll be able to perform examinations with greater ease and accuracy than ever before.

For more information on the Welch Allyn Elite Retinoscope, contact your Welch Allyn representative today or call our Customer Direct Center at 1-800-535-6663.

BRIGHTNESS VS. STANDARD RETINOSCOPE



Elite Optical System

The Elite Optical System provides 2.8X the brightness and superior streak quality versus standard retinoscopes.

US Patents

5,189,449 • 5,200,772 • 5,202,710

Lamp Patent

5,359,262



Spot Retinoscopy Interchangeability

Simply insert a spot filament bulb in place of the streak filament bulb, and the Welch Allyn Elite Retinoscope is converted into a Spot Retinoscope. This allows you to perform spot retinoscopy without having to buy a second scope.

Ordering Information

18245	Welch Allyn Elite Retinoscope
71000-C	Convertible Rechargeable Power Handle
71670	3.5v Rechargeable Handle for Desk Charger
71900	Lithium Ion Power Handle with AC Charging Module
71930/71932	Lithium Ion Handle with Ophthalmic Charging Stand Adapter
08200	Halogen HPX™ Streak Replacement Lamp
08300	Halogen HPX Spot Replacement Lamp

Superior Streak Quality

Enable easier and quicker refractions – even in patients with the smallest pupils

Crossed Linear Polarizing Filter

Eliminates 99% of bothersome glare from trial lenses, helping you perform a more precise exam at the flick of a switch

External Focusing Sleeve

The continuously rotating external focusing sleeve improves comfort and maneuverability of the scope, no matter what size your hand is, or how you prefer to manipulate the sleeve

Dust-Free Design

Special construction keeps the optics cleaner for more effective examinations



WelchAllyn®

Advancing Frontline Care™