



Introducing the New
Welch Allyn Vital Signs Monitor
300 Series

Advanced Monitoring Made Affordable



WelchAllyn

Advancing Frontline Care™

Welch Allyn Vital Signs Monitor 300 Series.

Simple. Flexible. Affordable.

Your spot-check and continuous vital signs monitor, all in one.

As a clinician you are under pressure to improve outcomes, boost operating efficiencies and lower costs every day. It's a daunting challenge, but one that can be made a little easier thanks to the new Welch Allyn Vital Signs Monitor 300 Series—an affordable vital signs monitor that's both flexible and easy to use.

Like the original Welch Allyn Vital Signs Monitor, the 300 Series helps enhance your productivity by providing automated vital signs before, during and after procedures. The simple, push-button menu operation and bright display allows you to quickly learn and easily obtain the important information you need, so you spend more time with your patients and less time with your equipment.

In addition to the great features you've come to expect from Welch Allyn, the Vital Signs Monitor 300 Series provides these exciting new features:

- Separate patient modes with individual settings for adult, paediatric and neonatal support
- Fast, accurate temperature readings with industry leading Welch Allyn SureTemp® Plus thermometry
- Simple, straightforward operation and intuitive menu navigation with our innovative monitor design

And because the 300 Series monitor is ideal for both routine assessment and continuous monitoring of vital signs, it's like getting two devices at one low price.



Procedures

The Vital Signs Monitor 300 Series is the ideal solution for taking routine patient vital signs or monitoring before during and after procedures, minor surgery and treatments in a variety of settings, including:

- Medical/Surgical Floors
- Paediatric Floors
- Emergency Departments
- Hospital Special Procedures Departments
- Labour and Delivery Floors
- Outpatient Surgery Centres
- Outpatient Clinics
- Specialty Physician Offices
- Urgent Care Centres
- Dialysis Centres
- Dental Offices

Comprehensive Features

Non-invasive Blood Pressure

The Welch Allyn Vital Signs Monitor 300 Series begins as a non-invasive blood pressure (NIBP) monitor. Even in its most basic configuration, it delivers the features you need.

- Multiple Patient Modes – Adult, Paediatric and Neonatal operating modes
- Automatic Blood Pressure Mode – choose from 13 measurement intervals
- Programmable Alarms – for high/low systolic, diastolic and pulse rate
- LCD Display – Clock window and menu navigation
- LED Display – Bright, clear viewing of vital sign readings
- Built-In Memory – save up to 99 sets of measurement data

4-Second SureTemp® Plus Thermometry

Add SureTemp Plus Electronic Thermometry and take advantage of the revolutionary measurement capability developed by Welch Allyn. Obtain oral, rectal and axillary temperatures with the speed you want and the accuracy you demand.

- Oral – approx. 4 seconds
- Axillary – Adults, approx. 15 seconds; paediatrics, approx. 10 seconds
- Rectal – approx. 10 seconds



Pulse Oximetry

The Vital Signs Monitor 300 Series can be configured with motion tolerant pulse oximetry (SpO₂) technology from Nellcor®. Benefit from accurate readings despite shaking or low perfusion.

- Programmable alarms – high/low saturation levels
- Wide variety of patient sensors
- Optional pulse tone feature – audible heart rate and SpO₂ saturation indicator

Optional Printer

Add an integrated Thermal Printer for trend analysis or record keeping and get what you need at the touch of a button.

Simple and Straightforward

The Welch Allyn Vital Signs Monitor has a very short learning curve. Simple push-button programming and a large, digital display make learning and using the monitor easy and hassle-free.

Choose Only the Features You Need

The Vital Signs Monitor 300 Series is designed so you can select precisely the features you need today. That means you don't have to pay for functions you don't use. If your monitoring needs change, upgrading is easy. The 300 Series supports upgrades to all functions – pulse oximetry, 4-second SureTemp Plus oral thermometry and integrated printer – so you can count on your monitor to accommodate your changing needs.





Take It Wherever You Need It

Carry it, set it on a countertop, wheel it on its own dedicated mobile stand, or mount it on a wall. Do it all—switching from one flexible mounting position to another is simple.

Vital Signs Monitor 300 Series goes wherever your patient goes—effortlessly.

Ordering Information:

53000	NIBP only (includes Pulse Rate and MAP)
5300P	NIBP/Printer
530T0	NIBP/Temperature
530TP	NIBP/Temperature/Printer
53N00	NIBP/Pulse Oximetry
53N0P	NIBP/Pulse Oximetry/Printer
53NT0	NIBP/Pulse Oximetry/Temperature
53NTP	NIBP/Pulse Oximetry/Temperature/Printer
008-0835-00	Mobile Stand with Accessory Basket
008-0834-00	Wall Mount with Accessory Basket

Language version suffix:

E=English, F=French, G=German, I=Italian, S=Spanish

Transformer version suffix:

1=North American, 2=European, 3=UK, 4=Australian

Welch Allyn

Advancing Frontline Care™

Welch Allyn UK Ltd, Cublington Road, Aston Abbots, Buckinghamshire HP22 4ND, England

Tel: + 44 (0) 207 365 6780 Fax: + 44 (0) 207 365 9694

Welch Allyn European Customer Services, Navan Business Park, Dublin Road, Navan, Co. Meath, Republic of Ireland

Tel: +353 46 90 67790 Fax: +353 46 90 67754

Welch Allyn SA, 394 Surrey Avenue, Ferndale 2194, P O Box 809, Fontainebleau 2032, South Africa

Tel: +27 11 777 7555 Fax: +27 11 777 7556