

## The Ultimate in Instrument Design

The new Pocket Professional Otoscope and Ophthalmoscope from Welch Allyn have been designed to provide you with instruments which are ideally suited for use both in the consulting room and when you are on the move. Extensive research was carried out during the development process to ensure that the size, shape and appearance of the instruments would meet the needs of today's busy health professionals - and at a competitive price.



Please complete your details on the reverse of this page and post or fax to your usual supplier for further information.

## Pocket Professional Diagnostic Set

The Pocket Professional Otoscope and Ophthalmoscope Diagnostic Set provides both of these new, modern instruments either in a presentation box or in a handy, durable carrying case.



Otoscope & Ophthalmoscope set:  
Code: 728359

Otoscope & Ophthalmoscope set in carrying case:  
Code: 728361

Part No. 03/0023-EUR



The New Pocket Professional  
Otoscope and Ophthalmoscope -  
The Ultimate in Instrument Design

**WelchAllyn®**

Advancing Frontline Care™

## Pocket Professional Otoscope

As you would expect from Welch Allyn, the new Pocket Professional Otoscope provides both the advantage of halogen light and a fibre optic light transmission system. This allows you to target cool, bright, white light exactly where you want it, leading to improved diagnosis and ease of use. But the Pocket Professional takes otoscope design into new territory - extensive research among physicians and other health professionals has been used to develop a beautifully designed handle which brings a new level of comfort, grip and feel - allowing it to sit well in the hand and be more convenient to use. Its sleek, modern design will also make an attractive addition to every consulting room.

- Ergonomically-designed handle sits comfortably in the hand and provides a secure grip
- Attractive, modern appearance
- Sealed system with insufflator port for pneumatic otoscopy
- Accepts disposable KleenSpec (or reusable SofSpec) specula
- Easy-to-use on/off switch

Available as single instrument, or in handy carrying case with space to add ophthalmoscope if desired. Takes 2 AA batteries (included, non-rechargeable). Comes with 4 disposable specula.

Otoscope:  
Code: 728351

Otoscope with carrying case:  
Code: 728353



## Pocket Professional Ophthalmoscope

Just like the Pocket Professional Otoscope, the new Pocket Professional Ophthalmoscope provides the advantage of halogen light coupled with precision optics, allowing you to target cool, bright, white light exactly where you want it, and providing the optimum conditions for accurate diagnosis and ease of use. Once again, Welch Allyn consulted regular ophthalmoscope users to develop an attractive, up-to-the-minute design which sits comfortably in the hand and allows that all-important sound grip, making it easy to use and manipulate.

- Ergonomically-designed handle sits comfortably in the hand and provides a secure grip
- Attractive, modern appearance
- Lens selection disc offers choice of 48 viewing-lens combinations
- Rubber brow rest prevents scratching of spectacles
- Sealed optics for improved long-term performance
- Continuous aperture selection dial offers six apertures (micro, small & large spot, fixation, slit & cobalt blue)
- Red free filter to ease identification of vessels and nerve fibres
- Cobalt blue filter can be used with fluorescein dye to quickly and easily detect corneal abrasions and foreign bodies
- Illuminated lens indicator makes lenses easier to select
- Easy-to-use on/off switch

Available as single instrument or in handy carrying case with space to add otoscope if desired. Takes 2 AA batteries (included, non-rechargeable).

Ophthalmoscope:  
Code: 728355

Ophthalmoscope with carrying case:  
Code: 728357



## Response Card

Please contact me regarding the Welch Allyn Pocket Professional Ophthalmoscope/Otoscope/Diagnostic Set (delete as applicable):

Name: .....

Address: .....

.....


..... Post Code: .....

Telephone: .....

Fax: .....

E-mail: .....

Please detach this portion and post/fax to your local distributor shown below:



Alternatively, please detach this portion and post/fax directly to Welch Allyn at the following address:

**Welch Allyn Limited**  
Navan Business Park, Dublin Road, Navan, Co Meath,  
Republic of Ireland

Tel. 00-353-469067790

Fax. 00-353-469067754

Welch Allyn (Europe) 2003

Part No. 03/0023-EUR

**WelchAllyn®**

Advancing Frontline Care™