

PURLEY

Fixed Height Couch



Strong and Stable Couch Frame with
Optional Drawer Units

www.sidhil.com



Features & Benefits

- Strong and stable couch frame construction with a fixed seat height of 760mm (30")
- Multi-position adjustable backrest for flexible patient positioning
- 50mm thick CMHR foam top for patient comfort
- Seamless rounded upholstery corners for hygiene and cleanliness
- Finished in durable magnolia epoxy coating with Biocote®
- Up to three optional drawer units can be fitted for storage of ancillary equipment
- Paper roll holder included which can be located at the head or foot end of the couch
- Quick and easy to assemble - tools supplied
- Fire retardant upholstery (BS5852) in seven colours
- One year warranty for peace of mind (excludes upholstery and product misuse)

UPHOLSTERY COLOURS



SPECIFICATIONS

Product Code	Description	Dimensions & Weights
COU31/Colour (*)	Purley Couch in epoxy magnolia	760mm(H) x 600mm(W) x 1830mm(L) Product weight: 31kg
COU96	Drawer Unit for Purley Couch (up to three can be fitted per couch)	136mm(H) x 562mm(W) x 593mm(L) Product weight: 14kg
COU71	Lithotomy Stirrups & Leg Holders (pair)	860mm(H) x 70mm(W) Product weight: 4kg
COU76	Replacement Holders for Stirrups	25mm(W) x 360mm(L) Product weight: 0.4kg
COU73/MAG	Single Tier Step in epoxy magnolia	155mm(H) x 530mm(W) x 310mm(L) Product weight: 3.4kg
COU78/MAG	Double Tier Step in epoxy magnolia	305mm(H) x 530mm(W) x 310mm(L) Product weight: 7.5kg

* Add name of colour to product code

All weights and measurements referred to in this leaflet are approximate. For exacting details contact the Sidhil Sales Department or ask your Sidhil Representative for further details.

The company reserves the right to change specifications without notice.

To see the full range of Doherty products offered by Sidhil, contact the Sales Support Team on
Tel: 01422 233000



Sidhil Business Park, Holmfild, Halifax
West Yorkshire HX2 9TN

Tel: +44 (0)1422 233000 Fax: +44 (0)1422 233010
Email: sales@sidhil.com www.sidhil.com



(93/42/EEC)



Certificate No. FM 14550



A member of the **Siddall & Hilton Ltd.** Group of Companies