

MERLIN MEDICAL

Warecrest C282 Rechargeable Cautery

DISCONNECT THE
MAINS SUPPLY BEFORE
REMOVING THIS COVER

PLEASE READ THE INSTRUCTION BOOKLET PROVIDED WITH THIS INSTRUMENT

CHARGE
OFF

MAINS
OFF

CE
0120

THE
CLASSIC
C282

Warecrest C282 Rechargeable Cautery

Warecrest



The new **C282** rechargeable cordless cautery is specifically designed to simplify minor surgery procedures. The C282 has many outstanding features & there is an extensive range of burners to choose from,

enabling the user to carry out a broad spectrum of procedures effectively & efficiently. The C282's portability & incredible power make it unique.

Established world wide the C282 is renowned by its users as powerful, practical & easy to use. Its uncompromising design means it can be situated anywhere in the room, erradicating the danger of hanging cables & wires.

The C282 is manufactured to meet the very highest demands of a wide range of users in science & medicine.

● Portability With Power

The C282 uses the very latest multicell technology (MCT) for its 2Ahr source. This means more power, giving more heat energy for longer periods of time. It is lightweight & portable & yet has the power of a mains operated unit.

● Simplicity With Comfort & Safety

With no trailing wires surrounding the user, treatment can be carried out safely in the most comfortable position for both doctor & patient & in almost any situation.

● Power At Your Fingertips

A simple operation requires just seconds of high heat energy power, yet the C282 can deliver several minutes of high temperature output, with just a few hours of recharging.

● Extensive Range Of Applications

The C282 can be used daily. The list below gives a general guide to some of its applications:

Skin tags, granuloma pyogenicum, subungual haematoma, accessory digits, entropion, warts, verrucae, condylomata, cervical erosions, spider naevi, squamous celled carcinoma etc.

The C282 is used by many other disciplines such as veterinary surgeons, chiropodists, laboratory technicians, animal welfare & industry.

● Cautery Ordering Information

Your C282 Rechargeable Cautery is supplied complete with battery charger, handle & rechargeable battery. A comprehensive user manual, protective cover & fitted mains plug (UK only) are also provided.

DESCRIPTION

220/240V AC 45/60Hz
110/120V AC 45/60Hz

ORDERING CODE

C282/220
C282/110

● Reliability At Its Highest

Because of its unique design & our attention to detail at the manufacturing stage the C282 is virtually maintenance free.

● The Modern Surgery

The C282 has been designed in such a way that it is practical yet stylish. Its compact size means it is always to hand & it will blend into its surroundings.

● Burners That Last

Every cautery needs burners in order to function. We use only precious metals to bring you the highest level of quality & reliability. In fact all our burners are handmade, ensuring a long life when used correctly.

● A Practical Range Of Burners

LENGTH		MODEL C282 & C28 BURNERS		BURNER SFTS													
mm	INches			A	B	C	D	E	F	G	H	J	K	L	M	X	Y
150	6"	RC00	CONCENTRATION LOOP														
120	5"	RC01	FLANGING POINT														
100	4"	RC02	RING COILS														
80	3"	RC03	STANDARD COILS														
60	2"	RC04	CO-AGULATION TAIL														
40	1"	RC05	COLD POINT														

Warecrest Website: www.warecrest.com

MERLIN MEDICAL

Distributed worldwide by:

**Merlin Medical Ltd.,
The Whitbread Centre,
Rhymney, Gwent, NP22 5XD, U.K.**

**Tel: +44 (0)1685 843 676
Fax: +44 (0)1685 843 860**

**E mail: sales@merlin-medical.co.uk
Website: www.merlin-medical.com**